

Kyaniusa Saalfori (Sterility Mist)

The Kyaniusa Saalfori an aerosol spray that is used by the [Poku Saeruo Degonjo](#) to decontaminate an area, it became available in 752 CY [YE 33](#).

About the Kyaniusa Saalfori

The Kyaniusa Saalfori is a potent aerosol vapor that kills any airborne or surface pathogen. Typically one or more cylinders are placed in a room and then activated with the vents shut down. Fifteen minutes of proper exposure is sufficient to ensure sterilization.

Warning

Personnel are advised to wear a filter mask or be in a sealed suit if they must be in the chamber. Breathing the vapor for a short duration can cause respiratory distress. Prolonged exposure will cause distress and may result in death. Prolonged exposure to skin may cause irritation, but is otherwise non-harmful.

General

Class: Sterilizing Type: Aerosol Designers: [Otâmovi Wiy Jael \(Silver Moon Sect\)](#) Manufacturer: [Otâmovi Wiy Jael \(Silver Moon Sect\)](#) on [Poku Movi Vamâqi \(Clan Void Berth\)](#)


Usage

Standard, is to place a cylinder in a room. One cylinder can fill a 5 meter cubed area. The nozzle has a built in delay timer. Turning the nozzle counter clockwise arms the cylinder. The user can select up to fifteen minute delay. Once the timer reaches zero the aerosol is released.

Typical decontamination uses are:

- Cargo hold after lifting from planet.
- Medical facility after contamination

Note: Kyaniusa Saalfori has recently be put to use in the [Anoka Saalfofi \(Sterility Shield\)](#)



From:
<https://wiki.stararmy.com/> - **STAR ARMY**

Permanent link:
https://wiki.stararmy.com/doku.php?id=faction:hidden_sun_clan:medical:sterility

Last update: **2023/12/21 04:22**

