Type 45 Standardized Field Medical Kit

Developed by Advancer Enterprises' Mechanics Division under the White Lion Brand, with consultation with the Yugumo Corporation, the Type 45 Standardized Field Medical Kit is designed to give medical practitioners the tools to operate in away from hospitals. It became available in YE 45.

Year of Creation	YE 45					
Designer	Yugumo Corporation, Advancer Enterprises					
Nomenclature	AE-G5-1A					
Alt. Version						
Manufacturer	Manufacturer Advancer Enterprises, Yugumo Corporation (Star Army version)					
Fielded by	Yugumo Corporation, Star Army of Yamatai and Customers					
Availability	Mass Production					
Price	800.00 KS 800 KS					

History

In YE 45, Yugumo Corporation, recognizing a market in need of up-to-date medical kits, approached Advancer Enterprises with a proposal to develop a comprehensive list of medical kits. One of them was a Standardized Field Medical Kit tailored for emergency trauma situations. Recognizing the opportunity for innovation and the potential to save lives on a grand scale, Advancer Enterprises' Mechanics Division eagerly accepted the challenge.

The project was initiated under Advancer's White Lion Brand, an up-and-coming subsidiary for highquality medical equipment. A cross-disciplinary team of engineers, medical experts, and field medics from both corporations was assembled to ensure the kit would meet the diverse needs of emergency personnel and the general public.

After months of rigorous research, prototyping, and field testing, the Type 45 Standardized Field Medical

Kit was finalized. The team faced numerous challenges during the development phase, including sourcing materials that would be both lightweight and durable, as well as ensuring the kit would be user-friendly for medical professionals and laypeople alike.

Kit Description

Beyond a simple trauma kit, the Type 45 Standardized Field Medical Kit is designed to allow medical personnel to perform a wide variety of medical functions in the field. While most of the contents lean towards trauma medicine the kit is closer to a doctor's bag for austere environments intended for use at least a terrain feature away from the front line.

The components are stored in a 60x40x30cm carrying case with a Kinugoshi-ko impregnated with Omnihue exterior and Yarvex lined interior with a thin layer of Yarvex Foam insulation sandwiched in between. Two large 30x15x10cm pounches can be attached to the side of the Standardized Field Medical Kit for expanded space.

The core kit can be accessed from the top compartment while splints and other durable items are held in the bottom compartment. Items can be either accessed from the top (protected with a drawstring and flap secured with buckles) to pull out individual component sections (after delatching it) or the front side of the kit can be opened as well, allowing the individual component sections to fold out.

Core	e Medical Components			
		PPE/Disinfectant/woun	d cleaning	
QTY	ltem	Trauma Type	Description/Use	Other Notes
20	Disposable Gloves	PPE	Disposable non-porous gloves with an antiseptic coating	
1	Sharps and Biohazard Containment Sleeve	Biohazard Cleanup	1.5L cylindrical sleeve designed for the safe disposal of sharps and medical waste	
1	Bottle of Antiseptic Skin Cleaner	Biohazard Cleanup/Sterilization	50mL Bottle of Antiseptic Skin Cleaner	
3	Biohazard Waste Bags	Biohazard Cleanup	Package of 3 bags for waste disposal.	
2	Polypropylene Biohazard Absorb Pads	Biohazard Cleanup	100L Absorption Value, a Polypropylene Pads 35cmx45cm folded.	
3	Surgical Masks	PPE ¹⁾	Individually Packaged Sterile a Surgical Masks	
1	Reusable Plastic Eye Shield	PPE ²⁾	A plastic shield to protect the eyes of the medic.	
1	Bottle of SEyewash	Clean Eyes and Eye Regions	Saline Eyewash solution.	
1	Bottle of Saline	Wound Irrigation	500mL bottle of Saline to clean wounds.	
4	Multipurpose Cloth	Biohazard Cleanup/Misc.	Cloths	

QTY			т	rauma Type		D	escrip	otion/Use		Other Notes	
10	Dispo Probe	sable Splin s	ter	Removal Objects			neld disposable probes, packages of 5 each.				
2	Emes	is Bag		Vomit Co	ontainment		1000mL Dis	posab	le Vomit Bag	js	
	-	c Gear									
QTY		ltem	Traum	а Туре		-	ion/Use		Othe	r Note	s
1	Traun Instru Bookl	ction	N/A		Additionally i to the latest medical data	KAIMO					
1	Biomonitor Bio			ics	Available wit KAIMON or P				Can also be KAIMON/Asc Companion		
20	Biopa	d	Biometr	ics	Cheap, adhe device that c and pulse				Good for pas observation	ssive	
1	ຈເຍ	penlight	Illumina		Pen-sized ligl illumination.	ht for lo	ocal area		Rechargeable Power Cells (5 pieces).		er Cells
10			and "Bit	Control e"	Depress tongue to allow mouth and throat examination.						
<u> </u>		s Prevent			sing Compo	1					
Q	YTY		I	tem		Trau	ита Туре		ription/Use	Othe	r Notes
	2	Gauze Eye	pads - C	val		Eye Coverings Eyer		Eyepa	ed Gauze ads, age of 2.		
	2				pressure	essure Excessive ble Bleeding an		used bleed and e	ned to be where ing is heavy xcessive. 12cm roll.		
	1				Wound Treatment Gel		90mL Medic				
	5	Advancer Enterprises Type 45 S Tourniquet			45 Smart	1	ire/Blood Control	Tourn	able a iquet easy ap-on/snap	Self-tig with d displa	0
	1	Clear 🔊 Su	urgical Ta	аре			ire-Sensitive ive Tape	-	n x 10m roll nesive tape.		
Multiple Rolls and Packages		and Bandages with Dispensers				Banda	ges	allows Banda Selec are In	tor Roll s for age Size tion. Some dividually aged, ³⁾		

PPE/Disinfectant/wound cleaning

Core Medical Components

		s Prevention/Wo			•	1	December			
	TY		Item			Trauma Type	Descriptio	n/Use	Other Notes	
Roll	ltiple s and kages	d Gauze Pads			Absorption		Various sizes of Gauze pads			
	Various sizes Adhesive, vented chest se			eals		Open Chest Wounds	One way aii	rflow	Allows air to escape the chest cavity but not blood	
5		S Abdominal/Combine Trauma Pad				An applied dressing to promote healing and protect the wound from further harm. Various Sizes.				
2 Molecu		Molecure Tape	Aolecure Tape			molecular bonding adhesive tape			Uses ranging from medical to utility	
		Absorbable sutures			X	X		X		
XAnti-clotting/ agent1Blutlassen		Anti-clotting/Antił agent	Antihemorrhagic/hemostati		atio	² X	x		x	
		Blutlassen			Х	X		X		
	1	Sinus Alleviation Spray			X	X		X		
	1	field surgical kit ir instruments	ncludes	Laparoscopic		X	delivery nai	nites	x	
	1	laparoscope				X X			X	
	1	small liquid-optim	ized va	cuum		X	Х		X	
Airw	/ays, l	Breathing, and C	irculat	ion Compor	ner	nts	•			
QTY		Item		Trauma Tyj	pe	Descriptio	n/Use	(Other Notes	
		nergency Pacemak	(er ⁴⁾	Absent heartbeat		Spike inserted into	sternum ⁵⁾	CPF	R alternative	
1	Auton	natic resuscitation	bag	Assessed breathing		A self contracting bladder Can be m pumped ⁶		n be manually mped ⁶⁾		
1	Small	Oxygen Tank		Х			(5-minute supply, can be connected to either of the above)		X	
1	Modular Oxygen Concentrator		Х		Х	X		X		
1	Compact Mechanical Ventilator		Х		Х	X				
X		Intubation tubes		Х		Endotracheal, reinfo			ko Sized	
Χ	Intubation tubes		Х		Endotracheal, reinforced, smooth Neple					
X		cal Cricothyroidoto	-	Х		Endotracheal, reinfo	orced, smoo	th Nej	pleslian Sized	
	1	l Pharma Compo								
QTY		Item		та Туре		Description/L			ther Notes	
1	╗Hya Gel	aluronic acid Burn	Burn			nL bottle of Burn Ge ling and hydration t	· ·	Inhibit Respo	s Inflammation	

2024/0	5/19 14.11
Burr	is and Pharma
QTY	Item

5/7

Burr	ns and Pharma Compo	onents		
QTY	Item	Trauma Type	Description/Use	Other Notes
5	Hyaluronic acid Burn Gel Packets	Burn	3.5mL Packet of Burn Gel, provides cooling and hydration to burns.	Inhibits Inflammation Response.
2	Burn Dressings	Burn	Silicone Gel lined Polyamide Burn Dressing. 10cmx10cm Individually wrapped.	
2	Burn Dressings	Burn	Silicone Gel lined Polyamide Burn Dressing. 5cmx15cm Individually wrapped.	
5	Hydrocortisone Anti- Itch Cream	rashes and skin irritation	3.5mL Packet of topical itch cream	
10	Acetaminophen	mild pain or fever	Pack of two 250mg chewable tablets	Do not take more than 1000mg at once or 4000mg in 24 hours ⁷⁾
10	lbuprofen	mild to moderate pain, inflammation, fever	Pack of two 200mg chewable tablets ⁸⁾	Do not take more than 800mg at once or 3200mg in 24 hours ⁹⁾
1	Nano-Enhanced Pain Relief Patches	Х	X	Х
Gen	eral Health Componer	nts		
QTY	Item		Trauma Type Description/Use O	ther Notes

QTY	Item	Trauma Type	Description/Use	Other Notes				
1	Portable Blood Analyzer	Х	X	Х				
1	Expandable Dental Emergency Kit	Х	X	Х				
1	Allergy and Anaphylaxis Kit	X	X	Х				
1	Personal Hygiene Pack	Х	Х	Х				
1	Environmental Adaptation Medication Kit	Х	Х	Х				
1	High Blood pressure medication	Х	Х	Х				
1	Glucose tablets	Х	X	Х				
1	Hypospray Injector	Х	X	Х				
2	Hypospray Cartridge, General Pain Relief	Х	Х	Х				
2	Hypospray Cartridge, General Pain Relief	Х	Х	Х				
1	Hypospray Cartridge, Anti-Inflammatory	Х	Х	Х				
1	Hypospray Cartridge, Internal Disinfectant	Х	Х	Х				
1	Hypospray Cartridge, Nolern	Х	Х	Х				
1	Hypospray Cartridge, Sedative	Х	Х	Х				
1	Hypospray Cartridge, Stimulant	Х	Х	Х				

Miscellaneous Gear

QTY	ltem	Trauma Type	Description/Use	Other Notes
4	Chemical Glowsticks	Local Illumination	A "crack" to activate chemical glow- light	
1	SAM Splint Roll	Splint Accessory	100cm roll of padded splint foam roll	

Miscella	neous Gear			
QTY	ltem	Trauma Type	Description/Use	Other Notes
1	1Surgical ForcepsRemoval of Medium-Small Foreign Objects		Sterilized Durandium Alloy forceps.	
1	Trauma Shears	Remove Clothing from Patients	Sterilized Durandium Alloy Shears.	
12	12 Safety Pins Clip and Hold		Titanium Alloy Safety Pins, 12 pc box.	
3	Oral Syringe	Liquid Delivery	12mL Oral Syringe, Disposable	
5	Butterfly Wound Closure Bandages	Bandages	Butterfly-style wound closure bandages, 5 individual packages.	
Various	JS Medical splint, and Alignment, and Placement Braces of Damaged Limb		Stored in the Bags bottom split case (removable)	
2	Instant Hot or Cold Compress	Pain or internal temperature regulation	1 x Instant Hot Compress, 1 x Instant Cold Compress. 10cm x 5cm packaged.	
5	5 Rescue Blanket Body Temperature Regulation		Mylar Rescue Blanket	
1	Compact Portable Shelter	Х	Х	Х

Type 45 Star Army Standardized Field Medical Kit

This is the Star Army of Yamatai version of the Type 45 Star Army Standardized Field Medical Kit. It includes all the same items as the Trauma Kit above, but has an product made by Ketsurui Zaibatsu and is considered restricted technology. These are carried by medical officers and can be refilled by Star Army Pharmacist.

	Type 45 Star Army	Standardize	d Field Medical Kit					
QTY	Item	Trauma Type	Description/Use	Other Notes				
Con	Contains all Items from the Type 45 Standardized Field Medical Kit Above Plus:							
IV Pouches	Includes a complete IV Kit	IV Kit ¹⁰⁾ and bags.	4 x 500mL bags of hemosynth ¹¹⁾¹²⁾					
2	Star Army Disposable Surgical Scrubs, Type 32	PPE	Cheap scrubs that can be disposed of to decrease cross contamination	Only in the Star Army Version				
1	Star Army Emergency Ration Pill		60 pill bottle					
1	Star Army Medication Anti- Radiation							
1	Star Army Body Bag, Type 35							
1	Star Army CASTER Tape Dispenser							

STAR ARMY - https://wiki.stararmy.com/

OOC Notes

2024/05/19 14:11

Locked Out created this article on 2022/12/06 08:36. Refined and submitted by Demibear here.

Products & Items Database				
Product Categories	kits, medical			
Product Name	Type 45 Standardized Field Medical Kit			
Nomenclature	AE-G5-1A			
Manufacturer	Advancer Enterprises, Yugumo Corporation			
Year Released	YE 45			
Price (KS)	800.00 KS			
Mass (kg)	20 kg			

¹⁾, ²⁾ Personal Protection Equipment

Multiple Types including Conforming Stretch Gauze Bandage, Compress Pressure Bandage, Triangular Bandage, and Elastic Wrap Bandage Roll

4)

Advancer Enterprises licensed version

5)

3)

Instructions included

Often called the external lung

7) 9)

if pain persists seek professional medical assistance

8)

take with food if possible

clamps, tubing, collars, spikes, drip chambers and vents

11)

This is for Nekovalkyrja only, non-differentiated hemosynth provides assistance and speeds up the natural Nekovalkyrja healing.

This product is not manufactured by Yugumo Corporation it is made by the Hemosynthetics division of Ketsurui Zaibatsu.



